

## Setup

Mode	ST
Flex control	yes
Bi-Flex	2
CPAP	12
AVAPS	yes
Vte	500
IPAP Max	25
IPAP Min	10
IPAP	20
EPAP	5
BPM	10
T <sub>i</sub>	2.0
Rise time control	yes
Rise time	3
Ramp time	10
Ramp start pressure	4
Tubing Type Control	on
Tubing Type	15
Apnea	10
Patient disconnect	60
Low Min vent	2.0
Low Vte	on
Backlight	on
Language	English
hPa/cmH <sub>2</sub> O	cmH <sub>2</sub> O
Setup Parameter Display	Vte
Humidifier	3
Reset Therapy Hours	no
Reset blower hrs	no
Provider mode	on
Back	