

SENSITIVITY FLOW PULSE

MODE

Control Manual Sigh

Control - Sigh

Assist Control

Assist/Control - Sigh

Assist

Spontaneous

SMV

Total Volume



.05 To 3 Liters

Rate Patient BPM



15 20 25

10 30 35

5 40

I:E Ratio



Blocks During Inverse I:E Ratio or Readings Over 1:9:9 (OFF During SMV)

Inspiratory Flow Rate



10 To 120 LPM

Increase Inspiratory Flow

PRESSURE

Sensitivity



Low Pressure Alarm Set

10 20 30 40

15 Sec. Delay

Exhalation Valve

Patient Air

PAP

Oxygen Sensor

Analyte Pressure Limit

OXYGEN



LIFECARE PLV-102

50

21

OFF

Percent

POWER

VAC

Internal Battery

External Battery

Battery Capacity

Reed Percent

All Flow Discontinuity

50% Recharge

Internal

External

Recharge

30 Second Alarm Silence

DANGER
RISK OF EXPLOSION IF USED NEAR FLAMMABLE ANESTHETICS